



Title: Inhaler Technique Review – Respimat Inhaler NMI-22

Patient det Title: Address:	Name: Re	ef number: rate of review:
GP, pharmacist, practice nurse or clinical specialist observes the patient completing the following steps and scores 1 (satisfactory) or 0 (not satisfactory) depending on whether the step is completed satisfactorily		
1	At time of dispensing, the pharmacist must ensure the patie demonstrate correctly loading and preparing the device or I prepare the device for first use for the patient. The device m discarded 3 months after first use.	load and
2	Hold the Respimat inhaler upright with the cap closed.	
3	Turn the base in the direction of the arrows until it clicks (ha	alf a turn).
4	Open the cap until it snaps fully open.	
5	In the sitting or standing position, with the head upright, bre gently (Do NOT breathe out into mouthpiece)	eathe out
6	Hold mouthpiece at right angles to the mouth and place bet teeth. Close lips around the mouthpiece and at start of insp (while inhaling slowly and deeply through the inhaler NOT to nose)	piration
7	press the dose release button and continue to breathe in a slong as possible. (Do NOT cover air vents with hands or m	
8	Remove the inhaler from the mouth and hold breath for 10 s as close to 10 seconds as possible	seconds, or
9	Breathe out slowly (Do NOT breathe out into the inhaler)	
10	Repeat steps 2-9 to receive a second dose, then close the ca	ар.
Total / 10		
	Signed/initialled by pharmacist/professional staff: Date of assessment:	
Recording of Patient Informed Consent: a) Patient's reason for choosing to attend for Inhaler Technique Review: b) I have received information about the Inhaler technique review process and give my consent: c) I agree that information may be shared with my GP or carer (specify) Date: d) I consent to the use of Inhaler Technique Review results, my name/identification having been removed, being used for audit purposes, including the accumulation of such data for research purposes.		
Patient signature:		Date:

Date: 01/9/15 Ver: 04